

NEW THERMAR

HEATING PAD



Targeted

The pillow is designed for a targeted application of heat. It is used for therapeutic purposes or simply as a source of immediate relief and relaxation for the body.



Adjustable

The control with display allows the selection of the desired temperature, with 4 levels of intensity adjustment.



Full body

It can be placed on multiple parts of the body (neck, shoulders, back, abdomen and joints).



Smart

Equipped with rapid heating and auto-off functions after 90 minutes.



Soft & igienic

Soft fabric lining washable at 30 °C.

INCLUDED ACCESSORIES

HEATING PILLOW

30 x 40 x 1 cm

POWER CABLE

Length: 2 m

CONTROL WITH DISPLAY

4 heat levels

PRODUCT CODE TM 100100

MAT-130

ANTI DECUBITUS MATTRESS



Adjustable bubbles type

Made of 130 air bubbles and equipped with end flaps for correct positioning in the bed. Pressure can be adjusted to best fit with patient weight.



Delicate

Massages skin delicately at alternating points of contact for improved blood flow and tissue oxygenation preventing bedsores formation.



Safe & non-toxic

Made of highly elastic non-toxic PVC, resistant and easy to clean.



Self inflating

Easily inflates with its own compressor.

INCLUDED PARTS

INFLATABLE MATTRESS

280 x 90 x 6,5 cm made of PVC

AIR PUMP

230V / 50Hz, Output: 5 l/min

PRODUCT CODE DE 620010